



POWERED BY
OLITEK




OzVader V1 Pandemic Ventilator


- Emergency use pandemic ventilator designed specifically to support global medical disaster situations
- Engineered to provide stable function in a wide variety of use locations (ie. medical disaster zones, hospital overflow treatment centers, during transport or regional medical clinics with limited hospital infrastructure)
- System does not require hospital pressured gas infrastructure or constant electrical power supply
- Flexible power systems including:
 - 110-240V AC
 - 12V DC
 - 12V internal battery
 - Air powered (medical or industrial)
- Intelligent robotic control system actuates a conventional Bag-Valve-Mask (BVM) to give precise control over key ventilation parameters.
- Automatically adapts to any adult BVM type (supplier agnostic)




OzVader – saving lives, anywhere.



- Medical Disaster Zones
- Hospital Overflow Treatment Centers



- Regional Medical Clinics



- Conventional Hospital Environments



COVID-19 Emergency Use Pandemic Ventilator

OzVader V1 Technical Specifications

Basic ventilation modes	Intermittent mandatory ventilation (IMV) Volume controlled ventilation (VCV)
Operational summary	<p>Robotic bag-valve-mask (BVM) actuation</p> <p>BVM actuated by the OzVader bellows system for long duration use</p> <p>Compatible with all adult sized BVM models (supplier agnostic).</p> <p>Set and monitor ventilation variables:</p> <ul style="list-style-type: none">▪ Tidal volume (Vt)▪ Respiratory rate▪ Flow rate▪ I:E ratio▪ PEEP▪ FiO2▪ Peak inspiratory pressure (PIP) <p>Backup internal air pump when compressed air is unavailable (12V)</p> <p>On-board micro-processor based controller with safety monitoring</p> <p>Simple user interface designed with ICU clinicians in mind for functionality & usability</p> <p>Flow and pressure sensors monitor ventilator and patient variables with integrated alarm functions</p>
Contamination control	<p>Removable patient exhalation module compatible with low temperature sterilization processes</p> <p>Compatible with viral filters on exhalation & inhalation limbs</p> <p>Compatible with invasive intubation</p>
Reliability & flexibility	<p>Designed for medical disaster situations requiring reliability, durability and flexibility with minimal moving parts and inbuilt power system redundancy</p> <p>Minimum 14 days continuous run-time on a single BVM</p>
Oxygen & FiO2	Recommended oxygen flow rate provided to user based on measured air flow. User sets the external wall/bottle mounted O2 flow control

	<p>between 0% and 100%</p> <p>Internal oxygen flow sensor monitors FiO2% with integrated oxygen visual and audible alarms</p>
Tidal volume	<p>Adjustable range of 200ml to 800ml (validated)</p> <p>Maximum range to 950ml (not formally validated)</p>
Respiratory rate	5 to 30 breaths per minute – adjustable
Flow	<p>5L/min to 100L/min – adjustable</p> <p>Up to 100L/min achievable when connected to external compressed air supply</p>
Inspiratory rate	0.5 to 2sec inspiratory time - adjustable
Inspiratory to expiratory time ratio (I:E)	4:1 to 1:9 ratio – adjustable
Peak inspiratory pressure (PIP)	<p>Direct measurement of peak inhalation pressure up to 80cm H2O - adjustable</p> <p>Mechanical fail-safe pressure valve at 80cmH2O maximum</p> <p>Target plateau pressures < 30cmH2O</p> <p>Driving pressure < 15cmH2O</p>
Positive end expiratory pressure (PEEP)	<p>Adjustable PEEP between 2-25cmH2O – adjustable</p> <p>Compatible with standard disposable & re-usable PEEP valves</p>
Ventilation circuit	<p>Twin tube ventilation circuits (inspiratory and expiratory limbs)</p> <p>Compatible with all standard 22 mm diameter (OD) 'male' standard circuit connectors</p>
ICU clinician functionality	<p>Clinical override function of 35cm H2O PIP cut-off while performing diagnostic procedures</p> <p>Inspiratory pause function allowing setting of safe tidal volumes and PEEP (plateau pressure measurement) – 10sec</p>
Alarmed functions	<ul style="list-style-type: none"> ▪ Power-off while ventilating ▪ External power supply failure ▪ Back-up battery ▪ Patient disconnect / patient circuit leak ▪ Low FiO2

- Maximum PIP exceedance
- Tidal volume not met or exceeded
- PEEP high/low

Power modes

110V-240V AC (plug type provided to suit delivery destination)

Internal 12V DC battery backup

External compressed air (hospital air / industrial compressed air)





Contact Details

Phone

OzVader

+61 7 3260 5390

Email

support@ozvader.com

Address

OzVader
248 Earnshaw Road
Northgate
Brisbane, Queensland, Australia

